

Product Name Report No. NT/MON/31-22/23 Extraction Ratio		P	RODUCT SPE	CIFICATION			
Product Name	ICNL/Botanical Name	NICOTIANA	A TABACUM	Date of Expiry		03 Years	
Report No.	Product Name	Nicotine Bit	artrate Dehydrate				
Parts of Plant Cas No 96055-45-7 SPECIFICATION METHOD TEST SPECIFICATION METHOD A. Physical Parameters: Vapearance dentification C by HPLC The retention time of major peak of chromatograph of the sample solution corresponds to that in the chromatogram of standard solution, as obtained in the assay test. Deganic impurities By HPLC As per USP Specification USP Vicotine Related Compound (A.B.C.D.E.F.G) NMT 0.3% USP USP NMT 0.1% USP NMT 0.1% USP NMT 0.5% USP NMT 8.0% w/w USP B. Active Ingredients: Vicotine Between 30.4% to 33.0% IPLC C. Heavy Metals: **ad (PB) **NMT 1.0 ppm **ICP-MS** Subminus (CD) **NMT 0.1 ppm **ICP-MS** Subminus (CD) **ICP-MS** Subminus (CD) **NMT 0.1 ppm **ICP-MS** Subminus (CD) **ICP-MS* Subminus (CD) **ICP	Report No.		<u> </u>				
TEST SPECIFICATION METHOD A. Physical Parameters: A. Physical Parameters: White to off white powder Visual The retention time of major peak of chromatograph of the sample solution corresponds to that in the chromatogram of standard solution, as obtained in the assay test. Deganic impurities By HPLC As per USP Specification USP Visual Organic impurities By HPLC As per USP Specification USP Victorine Related Compound (A.B.C.D.E.F.G) NMT 0.3% USP Vary Other Unspecified Impurities NMT 0.1% USP Cost on Drying (At 105°C) NMT 8.0% WW USP B. Active Ingredients: Victorine Between 30.4% to 33.0% HPLC C. Heavy Metals: Lead (PB) NMT 3.0 ppm ICP-MS Lorent (AS) NMT 1.0 ppm ICP-MS More	-					Indian NOM	
TEST SPECIFICATION METHOD A. Physical Parameters: A. Physical Parameters: A. Physical Parameters: A. Physical Parameters: White to off white powder dentification C by HPLC The retention time of major peak of chromatograph of the sample solution corresponds to that in the chromatogram of standard solution, as obtained in the assay test. Drganic impurities By HPLC As per USP Specification USP NNT 0.3% USP NNT 0.3% USP Osp Display (At 105°C) NNT 0.8% Osp Display (At 105°C) NNT 0.8% NNT 0.8% NNT 0.8% Second Pill Inpurities By HPLC NNT 0.8% NNT 0.8% NNT 0.8% NNT 0.8% NNT 0.8% ICP-MS Vivenic (AS) NNT 1.0 ppm ICP-MS Vivenic (AS) NNT 1.0 ppm ICP-MS NNT 0.1 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm NNT 0.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm NNT 1.0 ppm ICP-MS NNT 0.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 pcfu/g NNT 0.0 ppm ICP-MS NNT 0.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 pcfu/g NNT 0.0 ppm ICP-MS NNT 0.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate Count NNT 1.0 ppm ICP-MS D. Microbiological Test: Total Plate C							
A. Physical Parameters: Appearance White to off white powder Gentification C by HPLC The retention time of major peak of chromatograph of the sample solution corresponds to that in the chromatogram of standard solution, as obtained in the assay test. As per USP Specification USP Vicotine Related Compound (A,B,C,D,E,F,G) NMT 0.3% NMT 0.1% NMT 0.8% USP Loss on Drying (At 105°C) NMT 8.0% w/w USP B. Active Ingredients: NMT 0.8% NMT 8.0% w/w USP B. Active Ingredients: NMT 0.8 ppm ICP-MS Vicotine Between 30.4% to 33.0% HPLC C. Heavy Metals: Lead (PB) NMT 1.0 ppm ICP-MS ICP-MS ICP-MS IMMT 0.1 ppm ICP-MS ICP-MS ICP-MS IMMT 0.1 ppm ICP-MS ICP		90033-43-7	CDE CI				
Appearance White to off white powder White to off white powder Wrisual The retention time of major peak of chromatograph of the sample solution corresponds to that in the chromatogram of standard solution, as obtained in the assay test. Drganic impurities By HPLC As per USP Specification USP Nicotine Related Compound (A,B,C,D,E,F,G) NMT 0.3% USP NNT 0.1% NNT 0.1% USP Cost on Drying (At 105°C) NMT 8.0% w/w USP B. Active Ingredients: Nicotine Between 30.4% to 33.0% HPLC C. Heavy Metals: Lead (PB) NMT 3.0 ppm ICP-MS NNT 1.0 ppm ICP-MS Mercury (HG) NMT 1.0 ppm ICP-MS Mercury (HG) NMT 0.1 ppm ICP-MS D. Microbiological Test: D. Microbiological Test Cotal Plate Count NMT 10 cfu/g NMT 10 cfu/g Secherichia coli alaborat Absent USP 2022 Absent USP 2022 Secherichia coli alaborat Absent USP 2022 E. Mycotoxin: What is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material is completely required. During material handling, a full Lab PPE Kit should be kept away from sunshine & heat. The material is completely required. During material handling, a full Lab PPE Kit should be lested at all times. F. Other: PLEASE CARE FULLY READ THESE INSTRUCTIONS BEFORE OPENING THIS MATERIAL. The material is completely required. During material handling, a full Lab PPE Kit should be lested at all times. The material is product, Due to containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Bandard Packing: 25 Kg or 1 kg in single Container Customized Packing: 1 kg in single Container charges applicable Herball originar are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour			SPECI	IFICATION		METHOD	
The retention time of major peak of chromatograph of the sample solution corresponds to that in the chromatograpm of standard solution, as obtained in the assay test. Drganic impurities By HPLC As per USP Specification USP Note of the sample solution, as obtained in the assay test. Note of the sample solution, as obtained in the assay test. Note of the sample solution, as obtained in the assay test. Note of the sample solution, as obtained in the sample solution, as obtained in the assay test. Note of the sample solution of the sample solution in the sample solution of the sample solu	A. Physical Parameter	rs:					
of the sample solution corresponds to that in the chromatogram of standard solution, as obtained in the assay test. Drganic impurities By HPLC	Appearance		White to off white powder		Visual		
chromatogram of standard solution, as obtained in the assay test. As per USP Specification USP Nicotine Related Compound (A,B,C,D,E,F,G) NMT 0.3% NMT 0.1% USP Naty Other Unspecified Impurities NMT 0.1% USP Loss on Drying (At 105°C) NMT 8.0% w/w USP B. Active Ingredients: Nicotine Between 30.4% to 33.0% Between 30.4%	Identification C by HPLC		The retention time of major peak of chromatograph		USP		
the assay test. As per USP Specification USP Nicotine Related Compound (A,B,C,D,E,F,G) NMT 0.3% USP Nany Other Unspecified Impurities NMT 0.1% USP Cotal Impurities By HPLC NMT 0.8% USP NMT 0.8% USP Board Impurities By HPLC NMT 0.8% USP B. Active Ingredients: Vicotine Between 30.4% to 33.0% BPLC C. Heavy Metals:			of the sample solution corresponds to that in the				
Organic impurities By HPLC Nicotine Related Compound (A,B,C,D,E,F,G) NMT 0.3% USP NMT 0.1% IOSP NMT 0.1% IOSP NMT 0.8% USP B. Active Ingredients: NMT 0.8% NMT 8.0% w/w USP B. Active Ingredients: NMT 0.8% NMT 0.8% NMT 0.8% IPLC C. Heavy Metals: Lead (PB) NMT 3.0 ppm ICP-MS NMT 1.0 ppm ICP-MS ICP-MS			chromatogram of standard solution, as obtained in				
Nicotine Related Compound (A,B,C,D,E,F,G) NMT 0.3% USP NMY 0.1% USP NMT 0.1% USP Loss on Drying (At 105°C) NMT 8.0% W/W USP B. Active Ingredients: Nicotine Between 30.4% to 33.0% HPLC C. Heavy Metals: Lead (PB) NMT 1.0 ppm ICP-MS NAT 1.0 ppm ICP-MS Loss on Drying (At 105°C) NMT 1.0 ppm ICP-MS Loss on Drying (At 105°C) NMT 1.0 ppm ICP-MS Loss on Drying (At 105°C) NMT 1.0 ppm ICP-MS NAT 1.0 ppm ICP-MS Loss on Drying (At 105°C) NMT 1.0 ppm ICP-MS Loss on Drying (At 105			the assay test.				
Any Other Unspecified Impurities NMT 0.1% USP Fotal Impurities By HPLC NMT 0.8% USP B. Active Ingredients: NMT 0.8% W/W USP B. Active Ingredients: NOTE: NOTE	Organic impurities By HPLC				USP		
Cotal Impurities By HPLC	Nicotine Related Compound (A,B,C,D,E,F,G)				USP		
Cotal Impurities By HPLC	Any Other Unspecified Impurities				USP		
B. Active Ingredients: Nicotine Between 30.4% to 33.0% HPLC C. Heavy Metals: Lead (PB) NMT 3.0 ppm ICP-MS Cadmium (CD) NMT 1.0 ppm ICP-MS Cadmium (CD) NMT 1.0 ppm ICP-MS Cadmium (CD) NMT 0.1 ppm ICP-MS D. Microbiological Test: Fotal Plate Count NMT 10 cfu/g USP 2021 Sectorichia coli Absent USP 2022 Sectodomonas aeruginosa Absent USP 2022 Setudomonas aeruginosa Absent USP 2025 Setudomonas aeruginosa Absent USP 2026 Setudomonas aeruginosa Absent USP 2027 Setudomonas aeruginosa Absent USP 2028 Setudomonas aeruginosa Absent USP 2029 Setudomonas aeruginosa Absent USP 2020 Setudomonas 2020 Set	Total Impurities By HPLC						
Nicotine Between 30.4% to 33.0% HPLC C. Heavy Metals:	Loss on Drying (At 105°C)		NMT 8.0% w/w				
Nicotine Between 30.4% to 33.0% HPLC C. Heavy Metals:	B. Active Ingredients:				ı		
NMT 3.0 ppm ICP-MS	Nicotine		Between 30.4% to 33.0%		HPLC		
NMT 3.0 ppm ICP-MS	C. Heavy Metals:						
Arsenic (AS) NMT 1.0 ppm ICP-MS Addinum (CD) NMT 1.0 ppm ICP-MS ICP-MS ICP-MS Mercury (HG) NMT 0.1 ppm ICP-MS D. Microbiological Test: Total Plate Count NMT 100cfu/g Scherichia coli Absent USP 2021 Scherichia coli Absent USP 2022 Schamonella Absent USP 2022 Schapholococus aureus Et. Mycotoxin: Matoxin (B1+B2+G1+G2) Active Matoxin (B1+B2+G1+G2) The material is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot To avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour	Lead (PB)		NMT 3.0 ppm		ICP-M	S	
Cadmium (CD) NMT 1.0 ppm ICP-MS	Arsenic (AS)			1 1 1)) 1 1 1			
D. Microbiological Test: Fotal Plate Count NMT 100cfu/g USP 2021 Fotal Peace & Mold Count NMT 10 cfu/g USP 2021 Escherichia coli Absent USP 2022 Escherichia coli USP 2022 Escherichia (B1+B2+G1+G2) USP 561 Eschidual Pesticides USP 561 Eschidual Pesticides USP 561 F. Other: PLEASE CAREFULLY READ THESE INSTRUCTIONS BEFORE OPENING THIS MATERIAL. The material is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot fo avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Estandard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour	Cadmium (CD)			ALESIA V	ICP-M	S	
Total Plate Count NMT 100cfu/g USP 2021 Total Yeast & Mold Count NMT 10 cfu/g USP 2021 Escherichia coli Absent USP 2022 Salmonella Absent USP 2022 Seudomonas aeruginosa Absent USP 2022 Staphylococcus aureus USP 2022 Etaphylococcus USP 2022 Etaphyl	Mercury (HG)		NMT 0.1 ppm		ICP-M	S	
NMT 10 cfu/g	D. Microbiological Test						
Absent USP 2022 Seudomonas aeruginosa Absent USP 2022 Staphylococcus aureus Absent USP 2022 Staphylococcus aureus Absent USP 2022 Staphylococcus aureus Absent USP 2022 E. Mycotoxin: Matoxin (B1+B2+G1+G2) Assent USP 561 Residual Pesticides	Total Plate Count				USP 20	021	
Absent USP 2022 Staphylococcus aureus Absent USP 2022 Etaphylococcus aureus Matoxin (B1+B2+G1+G2)	Total Yeast & Mold Count		NMT 10 cfu/g				
Absent USP 2022 Staphylococcus aureus E. Mycotoxin: Matoxin (B1+B2+G1+G2) Assent USP 561 Residual Pesticides F. Other: PLEASE CAREFULLY READ THESE INSTRUCTIONS BEFORE OPENING THIS MATERIAL. The material is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot To avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour	Escherichia coli						
Absent USP 2022 E. Mycotoxin: Matoxin (B1+B2+G1+G2) USP 561 Residual Pesticides USP 561 F. Other: PLEASE CAREFULLY READ THESE INSTRUCTIONS BEFORE OPENING THIS MATERIAL. The material is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot To avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing:1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour	Salmonella						
E. Mycotoxin: Aflatoxin (B1+B2+G1+G2) USP 561 Residual Pesticides USP 561 F. Other: PLEASE CAREFULLY READ THESE INSTRUCTIONS BEFORE OPENING THIS MATERIAL. The material is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour							
Aflatoxin (B1+B2+G1+G2) Residual Pesticides USP 561 F. Other: PLEASE CAREFULLY READ THESE INSTRUCTIONS BEFORE OPENING THIS MATERIAL. The material is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot for avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour			Absent		USP 20)22	
Residual Pesticides USP 561 F. Other: PLEASE CAREFULLY READ THESE INSTRUCTIONS BEFORE OPENING THIS MATERIAL. The material is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot for avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour	•				,		
F. Other: PLEASE CAREFULLY READ THESE INSTRUCTIONS BEFORE OPENING THIS MATERIAL. The material is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot To avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour	,						
The material is completely required. During material handling, a full Lab PPE Kit should be used at all times. The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is ' Flushed ' and kept under a nitrogen ' Blanket ', and oxygen is evacuated from the bottle oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable. Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour				LICTIONS DEPONE O			
The material temperature should not exceed 25 °C throughout transit, processing, storage, should be kept away from sunshine & heat. The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot for avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable. Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour		A W ZIII A III				A W /III A I	
The substance ensures that if the bottle is opened, it is 'Flushed' and kept under a nitrogen 'Blanket', and oxygen is evacuated from the bot To avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour							
To avoid oxidation and degradation, the substance should be combined in an airless (Nitrogenous Environment). Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour							
Material mixing should take place in sealed containers. When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour							
When exposed to air, the material may quickly Oxidize; therefore, open the sample until it is ready to be tested or used, or develop your formulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable. Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour				i an anness (1410 genous Elli	, ii viiiil	ciitj.	
Cormulation. The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour				he sample until it is ready to 1	oe tested	d or used, or develop your	
The above substance meets with USP requirements. Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour	formulation.		,, open c			,	
Standard Packing: 25 Kg or 1 kg in single Container/Customized Packing: 1 kg in single Container charges applicable Herbal origins are used in the product. Due to seasonal and geographic variations in the plant material, there may be a slight Colour		USP requiremen	nts.				
	Standard Packing: 25 Kg or 1 kg	g in single Conta	iner/Customized Packin				
	Herbal origins are used in the pridifference.	roduct. Due to sea	asonal and geographic va	ariations in the plant material,	there n	nay be a slight Colour	

Note: This Is Computer Generated Specification & Does Not Require A Physical Signature.

